

Technical data sheet

Disposable Face Mask



Code:	ODZD-06043
Name of the product:	Disposable Face Mask
Producer:	For 2P SERVIS s.r.o., S PROMOTION s.r.o.

Description of the product:	For 2P SERVIS s.r.o., S PROMOTION s.r.o.
Material:	3 layers disposable face mask without sterilization
Size:	9,5 x 17,5 cm
Colour:	Blue, pink, green
Guideline: Standards: Certificates: Classification:	EN 14683:2019, Annex B, Bacterial Filtration Efficiency BFE \geq 98% ASTM F2101-19, Bacterial Filtration Efficiency BFE \geq 98% EN 14683:2019, Annex C, Differential Pressure to determine the breathability Delta P \leq 4 mm H ₂ O/cm ² ASTM F2101-19, Differential Pressure to determine the breathability Delta P \leq 40 Pa/cm ²
Using instruction:	<ul style="list-style-type: none"> Using instructions are not necessary because the correct use of the product is self-evident. Harmless if used by a properly qualified person. Use only undamaged products! Pay attention to the integrity of the product before and during use! The user is expected to test the suitability of the product for the intended purpose of use at his own risk!
Purpose and field of application:	It is suitable for health protection
Using restriction:	The product is disposable
General safety instructions:	After disposal, safe disposal is necessary to avoid environmental pollution It is not allowed to use in case of damaged inner packaging
Storage:	The product should be stored in a cool, dry, wellventilated and clean place with relative humidity of not more than 80%, without direct sunligh
Minimum durability:	Two years if stored under required conditions
Allergenicity information:	Do not use in case of allergy or sensitivity to any material of which the product is composed

Packing:	Amount:	
Packaging:	10 pcs	PE bag
	50 pcs	Box
	1000 pcs	Carton
Maintenance	The product is for single use, there is a risk of contamination with repeated use and deterioration of its properties	
Disposal:	The product can be thermally recycled or disposed of in landfills without releasing toxic substances. In case of contaminated products, disposal is determined by the nature and extent of the contamination in compliance with applicable laws and regulations	